

2

Repeat wt-1 peptide binding Assay.

5x10⁵ T2 cells/well in RPMI + 5% Boiled FCS.

Add peptides to a final conc. of 100 nM.

Leave to binding at 37°C 2h

wash x1 with PBS + 2% FCS.

Apply mAb A2 mAb HB54. HB117 30' on ice

wash x3

Apply G2M FITC (1/250) 30' on ice

wash x3

Resuspend cells in PBS (400ul).

FACS Analysis

Results, Good.

peptide

2nd Ab only

Mean: HB54

HB117

513

37.70

42.32

HPV E7/86-93

109.70

118.12

HRV

146.86

220.67

WT-1 /10-

~~93.04~~

144.38 #

WT-1 /17-

43.79

54.04 #

WT-1 /126-

69.78

30.79 #

WT-1 /187-

93.01

136.27 #

WT-1 /225-

70.03

97.56 #

WT-1 /235-

95.63

105.90 #

WT-1 /280-

45.66

76.42 +

WT-1 /441-

137.32

172.31 #

BEST AVAILABLE COPY

A2 Binding Titration

2

peptides

WT:

10

17

126

187

225

235

280

441

HPVE7 EG-93

HBV 17-28

Plate peptide in RPMI + 5% fetal FCS in 96 well plate (U bottomed), from 100 μ M of first column to column 11 (0.0001 μ M) in 50 μ l medium.

Add 50 μ l of T2 cells at a conc of 6×10^6 / ml i.e. 5×10^5 / well mix well.

37°C overnight.

wash x 1 \times PBC (cold) + 1% FCS

stain \times HB 117 and HB 54 (parent culture sp) 10 min for 30' on ice.

wash \times 3

Apply Ga M Ig FITC at 1/250 (50 μ l) for 30' on ice


wash \times 3

Analyse on FACS


Results: Good see attached sheets.

3

WT-1 peptides binding to A2 molecules on T2 cells
--- Stained with HB54 mAb



Binding of WT-1 peptides to A2 molecules on T2 cells
--- Stained with HB117 mAb



2 no pl

FL1



2 no pl



FL1

pop. HB



FL1-H

3.)

4

Setting up Allo- & Self restricted peptide CTL lines

MHC: L63 A2⁻ Allo

peptides: L64 A2⁺ Self

WT1 10, 235, 441 and C228 (true Control).

Method: as left use DC from L64 as stimulator

Stimulate with peptide loaded CIRAZ cells (50 μM).

Restimulated E " RMA-S A2 (control)

CTL assay

Expand & cloning

set up clones

Transfer to a well to forward plate

feed

5

Setting up: Self & Allo - restrict peptide specific CTL line.

PBMC	L G 5	A2+	Self
	L G 6	A2-	Allo

pep tiles: WT1 126, 187, 280
+ve control E786-95.

Methods, see previous page.

First round stimulation with T2 loaded peptides.

Restimulate with peptide loaded - CIR A2

4/2	cl	RMA-S - A2
5/2	"	"
7/3	Restimulation	CIR A2
7/3	CTL	
23/3	Cloning	
10/4	CTL of clones	

6

Protocol for peptide specific allorestricted human CTL

1. Prepare stimulator cells: in case of T2 stimulators preincubate o/n 10^6 cells/ml/well in x24 well plates in RPMI medium plus 5% boiled (10 min) FCS with 100 nM peptide. If allo DC are used as stimulators, allo DC from A2 +ve individuals can be separated as following:

Incubate A2+ PBMC at 37°C O/N, harvest nonadherent cells, resuspend in RPMI+9%FCS at 5×10^6 /ml and layer onto metrizadime gradient carefully, spin at 1500 for 10 min and harvest the interface (30% low density DC, 20% recovered from PBMC counts). Wash X3 with medium, resuspend at 10^6 /ml, add peptides to a final conc of 50uM and incubate at 37°C for 2 hrs. *3000 rad*
The stimulator can be added to responders directly w/O washing off peptides.

2. Separate blood by Ficoll centrifugation, wash twice, count, make sure that cells are HLA-A2 negative (by staining with HB54 and HB117 antibodies)

3. Plate 2×10^6 effectors and 2×10^5 stimulators per well in x24 well plates in culture medium (RPMI 1640, 10% FCS, 1% glutamine, 1% pen/stropt 5×10^{-5} M 2-ME) with 500 nM peptide

4. On day 4 prepare stimulators as described in 1. if required. It may be not necessary to preincubate C1RA2 cells for more than 2 hrs

5. On day 5 harvest cells, count spin and seed 5×10^5 effectors, 2×10^6 irradiated (about 3000 rad) HLA-A2 negative autologous feeders, 2×10^5 peptide loaded and irradiated (about 10000 rad) T2 or C1RA2 cells and 500 nM peptide in 2ml of the culture medium with 10u/ml IL-2, 2.5ng/ml IL7 (1/2000 of the stock of 5ug/ml) and with 10% of Q12054 culture sup (anti-CD4 antibodies) (complete medium)


6. 14 days (and further every 2 weeks) later restimulate cultures as described in 5. Cell recovery is usually 200%

7. On day 5 analyse peptide specificity in a ^{51}Cr release assay using T2 or mouse cells transfected with A2 molecule as CTL targets and cold K562 cells as NK targets (10 per 1 labelled target). Use 0.7% TFA for max release instead of 0.5% SDS. Significant peptide specificity is detectable after 5-7 week in culture with E:T ratio 10:1

8. Peptide specific cultures can be cloned seeding 1, 10 and 100 cells/well in round bottom x96 well plates in complete culture medium with 10^4





4

Sheet1



8

Sheet1



572	-9
	-16

67

19

test done 31 Co Release Assay
Standard 4 hr


target T2 + E.7
T2 + 126
K562
K562 Az

more tests


Sheet 1

10


hWT126 clone77




hWT126 clone 81




HWT126 CLONE85




HWT126 CLONE62



HWT126 CLONE 65



HWT126 CLONE1



771	-3	1199.2	4
907.4	0	979.2	0
	-9		74
	-2		79
	-2		25
	0		81

7)

11

51 Cr assay

clones 77. 81. 32.

targets: T2 + E7 vs. control

T2 + 126

Bv 173

UV 441

697

CMC Bone Marrow 4 (Ferguson)

CD34+ from BM (by Steve)

E = T 25. 12. 6. 3. 1.0 075


Results: line 77 line 32.

	line 77		line 32
		*	81
T2 + E7	+	+	+
T2 + 126	+	+	+
Bv 173	+	+	+
697	+	+	+
UV 441	+	+	+
CMC CD34+	+	+	-
CMC BM 4	-	+	-
CMC BM 4			

The clone 77 kills CD34+ from CMC BM

12

HWT12 clone 77 @CD34 CML BM



41
12
10
23
27
76
69
29
80
0
11
60
68
40
31
97
97
100
105
19
117
99
100
108
25
85
102
61
88
99
87
12
61
25
101
89
24
81
85
26
118
26
28
46
96
69
1
76

			-1	
			-4	
			-1	
			-1	
			28	
			57	
697	RNA	RT PCR	19	
BV173		+	50	
		##	-2	
Ferguson	cd34	+	-2	
Thomas		+	-4	
Grady		+	-2	
			-1	
			67	
			13	
			6	
			3	
			0	
			-5	
			-5	
			-7	
			13	
			-3	
			-5	

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

BLACK BORDERS

IMAGE CUT OFF AT TOP, BOTTOM OR SIDES

FADED TEXT OR DRAWING

BLURRED OR ILLEGIBLE TEXT OR DRAWING

SKEWED/SLANTED IMAGES

COLOR OR BLACK AND WHITE PHOTOGRAPHS

GRAY SCALE DOCUMENTS

LINES OR MARKS ON ORIGINAL DOCUMENT

REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY

OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.